

Anti-Cr2/Cd21 Rabbit Polyclonal Antibody For Wb, Ihc, If, Elisa - P20023

Item Number: CM0028269



Introduction

Anti-CR2/CD21 rabbit polyclonal antibody (Gene ID: 1380) for WB, IF, IHC, ELISA applications. Validated in human, mouse, rat. Recombinant immunogen targeting complement receptor type 2 (CD21), the Epstein-Barr virus receptor and B-cell coreceptor.

[Learn More](#)

Parameter	Value
Product Name	CR2/CD21 Rabbit pAb
Remarks/Alias	CR; C3DR; CD21; CVID7; SLEB9; CR2/CD21
Species	Human
GeneID (Human)	1380
GeneID	1380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-590 of human CR2/CD21 (NP_001006659.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20023
Protein Weight	113kDa
Shipping	Ice bag