

Anti-Crabp1 Monoclonal Antibody For Wb, If-P, Elisa - P29762

Item Number: CM0011409



Introduction

CamelBio offers anti-CRABP1 rabbit monoclonal antibody validated for WB, IF-P, and ELISA. Detects human, mouse, and rat CRABP1, a cytoplasmic retinoic acid-binding protein involved in retinoid signaling. Ideal for cellular retinoic acid studies.

[Learn More](#)

Parameter	Value
Product Name	CRABP1 Rabbit mAb
Remarks/Alias	RBP5; CRABP; CRABPI; CRABP-I; CRABP1
Species	Human
GeneID (Human)	1381
GeneID	1381
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 58-137 of human CRABP1 (P29762).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29762
Protein Weight	16kDa
Shipping	Ice bag