


Anti-Crabp2 Polyclonal Antibody For Wb, If/Icc, Elisa - P29373

Item Number: CM0027351



Introduction

CRABP2 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA applications. Detects CRABP2 in human, mouse, and rat samples. Ideal for retinoic acid signaling research.

[Learn More](#)

Parameter	Value
Product Name	CRABP2 Rabbit pAb
Remarks/Alias	RBP6; CRABP-II; CRABP2
Species	Human
Gene ID (Human)	1382
Gene ID	1382
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 22-94 of human CRABP2 (NP_001869.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29373
Protein Weight	16kDa
Shipping	Ice bag