

Anti-Cracr2A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bsw2

Item Number: CM0017120



Introduction

Rabbit polyclonal antibody against CRACR2A (RAB46/EFCAB4B) for WB and ELISA. Reacts with human, mouse, rat. Involved in store-operated Ca²⁺ entry, CRAC channel regulation, and T-cell signaling.

[Learn More](#)

Parameter	Value
Product Name	CRACR2A Rabbit pAb
Remarks/Alias	RAB46; EFCAB4B; CRACR2A
Species	Human
GeneID (Human)	84766
GeneID	84766
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CRACR2A (NP_001138430.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BSW2
Protein Weight	83kDa
Shipping	Ice bag