

Anti-Cracr2A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bsw2

Item Number: CM0014472



Introduction

High-quality rabbit polyclonal antibody against human CRACR2A, validated for WB, IHC-P, and ELISA. Detects endogenous CRACR2A in human and mouse samples, essential for studying store-operated calcium entry and T-cell activation.

[Learn More](#)

Parameter	Value
Product Name	CRACR2A Rabbit pAb
Remarks/Alias	RAB46; EFCAB4B; CRACR2A
Species	Human
GeneID (Human)	84766
GeneID	84766
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CRACR2A (NP_001138430.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BSW2
Protein Weight	83kDa
Shipping	Ice bag