

Anti-Cralbp Rabbit Polyclonal Antibody For Wb, Elisa - P12271

Item Number: CM0029831



Introduction

Rabbit polyclonal CRALBP antibody validated for WB and ELISA, reactive with human, mouse, rat. CRALBP is a retinoid carrier essential for visual cycle regeneration and vision research.

[Learn More](#)

Parameter	Value
Product Name	CRALBP Rabbit pAb
Remarks/Alias	CRALBP
Species	Human
GeneID (human)	6017
GeneID	6017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-317 of human CRALBP (NP_000317.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P12271
Protein Weight	36kDa
Shipping	Ice bag