

Anti-Crat Polyclonal Antibody For Wb, If/Icc, Elisa - P43155

Item Number: CM0030503



Introduction

Explore our ready-to-use CRAT Rabbit Polyclonal Antibody for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat CRAT (71 kDa) for metabolism and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	CRAT Rabbit pAb
Remarks/Alias	CAT; CAT1; NBIA8; CRAT
Species	Human
GeneID (Human)	1384
GeneID	1384
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-626 of human CRAT (NP_000746.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43155
Protein Weight	71kDa
Shipping	Ice bag