

Anti-Crb1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P82279

Item Number: CM0013457



Introduction

CRB1 rabbit polyclonal antibody validated for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Confirmed reactivity with human and mouse tissues. Ideal for retinal development and cell polarity research. SWISS accession P82279.

[Learn More](#)

Parameter	Value
Product Name	CRB1 Rabbit pAb
Remarks/Alias	LCA8; RP12; CRB1-A; CRB1-B; CRB1-C; CRB1
Species	Human
GeneID (Human)	23418
GeneID	23418
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-133 of human CRB1 (NP_957705.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P82279
Protein Weight	154kDa
Shipping	Ice bag