

Anti-Crbn Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Sw2

Item Number: CM0006415



Introduction

Rabbit monoclonal anti-CRBN antibody validated for WB, IHC-P, IF/ICC, ELISA. Targets cereblon, a substrate receptor in the E3 ubiquitin ligase complex, essential for proteasomal degradation. Reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	CRBN Rabbit mAb
Remarks/Alias	MRT2; MRT2A; CRBN
Species	Human
Gene ID (Human)	51185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CRBN (NP_057386.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96SW2
Protein Weight	51kDa
Shipping	Ice bag