


Anti-Crcp Monoclonal Antibody For Wb And Elisa - 075575

Item Number: CM0007412



Introduction

Rabbit monoclonal antibody targeting human CRCP (DNA-directed RNA polymerase III subunit RPC9 / CGRP-receptor component protein). Suitable for WB and ELISA; cross-reacts with mouse. Ideal for transcription, innate immunity, and CGRP signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CRCP Rabbit mAb
Remarks/Alias	C17; RCP; RCP9; RPC9; POLR3I; POLR3J; CGRPRCP; CGRP-RCP; CRCP
Species	Human
GeneID (human)	27297
GeneID	27297
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CRCP (NP_001035737.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	075575
Protein Weight	17kDa
Shipping	Ice bag