


Anti-Crcp Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 075575

Item Number: CM0027811



Introduction

Rabbit polyclonal antibody against CRCP (CGRP-RCP, POLR3I) for Western blot, IHC, and ELISA. Validated in human, mouse, and rat samples. Suitable for RNA polymerase III and CGRP receptor studies.

[Learn More](#)

Parameter	Value
Product Name	CRCP Rabbit pAb
Remarks/Alias	C17, RCP, RCP9, RPC9, POLR3I, POLR3J, CGRPRCP, CGRP-RCP, CRCP
Species	Human
Gene ID	27297
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human CRCP (NP_001035737.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075575
Protein Weight	17kDa
Shipping	Ice bag