

Anti-Creb1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P16220

Item Number: CM0015460



Introduction

Rabbit polyclonal antibody against human CREB1, validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Ideal for transcription regulation, circadian rhythm, and cell differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	CREB1 Rabbit pAb
Remarks/Alias	CREB; CREB-1; CREB1
Species	Human
Gene ID (Human)	1385
Gene ID	1385
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-341 of human CREB1 (NP_604391.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16220
Protein Weight	35kDa
Shipping	Ice bag