

Anti-Creb1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P16220

Item Number: CM0016579



Introduction

High-quality rabbit polyclonal antibody against human CREB1, validated in WB, IP, and ELISA. Useful for studying transcriptional regulation, circadian rhythms, and cellular differentiation. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	CREB1 Rabbit pAb
Remarks/Alias	CREB; CREB-1; CREB1
Species	Human
GeneID (Human)	1385
GeneID	1385
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human CREB1 (NP_604391.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16220
Protein Weight	35kDa
Shipping	Ice bag