

Anti-Creb3L2 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q70Sy1

Item Number: CM0005660

Introduction

Rabbit monoclonal antibody against human CREB3L2 (BBF2H7). Validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat CREB3L2 (57 kDa). Suitable for ER stress and chondrogenesis research.

[Learn More](#)

Parameter	Value
Product Name	CREB3L2 Rabbit mAb
Remarks/Alias	BBF2H7
Species	Human
GenelD (Human)	64764
GenelD	64764
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-247 of human CREB3L2 (NP_001240704.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q70SY1
Protein Weight	57 kDa
Shipping	Ice bag