

Anti-Creb3L2 Rabbit Polyclonal Antibody For Wb, Elisa - Q70Sy1

Item Number: CM0018293

Introduction

CREB3L2 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human and mouse CREB3L2, a 57kDa ER stress-responsive transcription factor involved in chondrogenesis and UPR.

[Learn More](#)

Parameter	Value
Product Name	CREB3L2 Rabbit pAb
Remarks/Alias	BBF2H7; CREB3L2
Species	Human
Gene ID (Human)	64764
Gene ID	64764
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 1-247 of human CREB3L2 (NP_001240704.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q70SY1
Protein Weight	57kDa
Shipping	Ice bag