

Anti-Creb5 Rabbit Polyclonal Antibody For Wb And Elisa - Q02930

Item Number: CM0013535

Introduction

CREB5 rabbit polyclonal antibody validated for WB and ELISA, suitable for human and rat samples. Targets transcription factor CREB5, which binds cAMP response element to regulate gene expression.

[Learn More](#)

Parameter	Value
Product Name	CREB5 Rabbit pAb
Remarks/Alias	CREB-5; CREBPA; CRE-BPA; CREB5
Species	Human
GeneID (human)	9586
GeneID	9586
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-380 of human CREB5 (NP_878902.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q02930
Protein Weight	57kDa
Shipping	Ice bag