

Anti-Crebzf Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ns37

Item Number: CM0028318

Introduction

CREBZF rabbit polyclonal antibody validated for Western blotting and ELISA. Detects human, mouse, rat CREB/ATF bZIP transcription factor, involved in transcriptional regulation and HSV suppression. Ideal for gene expression studies.

[Learn More](#)

Parameter	Value
Product Name	CREBZF Rabbit pAb
Remarks/Alias	ZF; SMILE; CREBZF
Species	Human
GenelD (Human)	58487
GenelD	58487
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NS37
Protein Weight	37kDa
Shipping	Ice bag