

Anti-Creg1 Rabbit Polyclonal Antibody For Wb, Elisa - 075629

Item Number: CM0024502



Introduction

Rabbit polyclonal antibody targeting human CREG1, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Recognizes a ~24 kDa protein; immunogen corresponds to amino acids 110-190. Ideal for studying transcriptional control and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	CREG1 Rabbit pAb
Remarks/Alias	CREG; CREG1
Species	Human
Gene ID (Human)	8804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-190 of human CREG1 (NP_003842.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075629
Protein Weight	24kDa
Shipping	Ice bag