

Anti-Creld1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Hd1

Item Number: CM0012892

Introduction

High-quality Anti-CRELD1 polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat CRELD1 (45 kDa), a protein disulfide isomerase involved in AChR localization. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	CRELD1 Rabbit pAb
Remarks/Alias	AVSD2; CIRPIN; CRELD1
Species	Human
GeneID (Human)	78987
GeneID	78987
Immunogen	Synthetic Peptide (a synthetic peptide corresponding to a sequence within amino acids 50-150 of human CRELD1 (NP_056328.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96HD1
Protein Weight	45kDa
Shipping	Ice bag