


Anti-Crh Rabbit Polyclonal Antibody For Wb And Elisa - P06850

Item Number: CM0018425



Introduction

Rabbit polyclonal antibody against human Corticotropin-releasing hormone (CRH), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for stress and neuroendocrine research.

[Learn More](#)

Parameter	Description
Product Name	CRH Rabbit pAb
Remarks/Alias	CRF; CRH1; CRH
Species	Human
GeneID (human)	1392
GeneID	1392
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-194 of human CRH (NP_000747.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06850
Protein Weight	21kDa
Shipping	Ice bag