

Anti-Crhbp Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P24387

Item Number: CM0025299

Introduction

Rabbit polyclonal antibody against human CRHBP, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for studying CRH/urocortin signaling and HPA axis regulation.

[Learn More](#)

Parameter	Value
Product Name	CRHBP Rabbit pAb
Remarks/Alias	CRFBP; CRF-BP; CRHBP
Species	Human
GeneID (Human)	1393
GeneID	1393
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-322 of human CRHBP (NP_001873.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24387
Protein Weight	36kDa
Shipping	Ice bag