


Anti-Crhr1 Rabbit Polyclonal Antibody For Wb And Elisa - P34998

Item Number: CM0028734



Introduction

Anti-CRHR1 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, rat CRHR1. Ideal for stress response, neuroendocrine signaling, and GPCR research. Uniprot: P34998.

[Learn More](#)

Parameter	Value
Product Name	CRHR1 Rabbit pAb
Remarks/Alias	CRF1; CRHR; CRF-R; CRFR1; CRF-R1; CRFR-1; CRH-R1; CRHR1L; CRF-R-1; CRH-R-1; CRHR1
Species	Human
GeneID (Human)	1394
GeneID	1394
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-120 of human CRHR1 (NP_001138618.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34998
Protein Weight	51kDa
Shipping	Ice bag