


Anti-Crip2 Rabbit Polyclonal Antibody For Wb, Elisa - P52943

Item Number: CM0025469



Introduction

Rabbit polyclonal antibody against human CRIP2 (ESP1), validated for Western blot and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying cysteine-rich protein 2 function in cardiac and other tissues.

[Learn More](#)

Parameter	Value
Product Name	CRIP2 Rabbit pAb
Remarks/Alias	CRIP; CRP2; ESP1; CRIP2
Species	Human
GeneID (Human)	1397
GeneID	1397
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human CRIP2 (NP_001303.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52943
Protein Weight	22kDa
Shipping	Ice bag