


Anti-Crisp2 Rabbit Polyclonal Antibody For Wb, Elisa - P16562

Item Number: CM0027535



Introduction

High-quality rabbit polyclonal antibody targeting human CRISP2 (CT36/TPX1). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for sperm capacitation and reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	CRISP2 Rabbit pAb
Remarks/Alias	CT36; TPX1; TSP1; GAPDL5; CRISP-2; CRISP2
Species	Human
GeneID (Human)	7180
GeneID	7180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-243 of human CRISP2 (NP_003287.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P16562
Protein Weight	27kDa
Shipping	Ice bag