

Anti-Crispld2 Polyclonal Antibody For Wb And Elisa - Q9H0B8

Item Number: CM0016112

Introduction

Rabbit polyclonal antibody against human CRISPLD2, validated for Western blotting (WB) and ELISA. Cross-reacts with mouse and rat. Suitable for studying CRISPLD2 function in matrix assembly and extracellular matrix biology.

[Learn More](#)

Parameter	Value
Product Name	CRISPLD2 Rabbit pAb
Remarks/Alias	LGL1; CRISP11; LCRISP2; CRISPLD2
Species	Human
Gene ID (Human)	83716
Gene ID	83716
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 140-430 of human CRISPLD2 (NP_113664.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H0B8
Protein Weight	56kDa
Shipping	Ice bag