


Anti-Crk/Crkii Polyclonal Antibody For Wb, Elisa - P46108

Item Number: CM0024625



Introduction

CamelBio offers anti-CRK/CRKII rabbit polyclonal antibody for Western blot and ELISA applications. It detects human and mouse CRK, suitable for studying cell adhesion, migration, and phagocytosis signaling. SWISS: P46108.

[Learn More](#)

Parameter	Value
Product Name	CRK/CRKII Rabbit pAb
Remarks/Alias	p38; CRKII; CRK/CRKII
Species	Human
Gene ID (Human)	1398
Gene ID	1398
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-304 of human CRK/CRKII (NP_058431.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46108
Protein Weight	34kDa
Shipping	Ice bag