


Anti-Crk/Crkii Rabbit Monoclonal Antibody For Wb, Elisa - P46108

Item Number: CM0009921



Introduction

High-quality rabbit monoclonal antibody targeting human CRK/CRKII (p38) for WB and ELISA. Validated in human and mouse samples. Ideal for cell adhesion, migration, and phagocytosis research.

[Learn More](#)

Parameter	Value
Product Name	CRK/CRKII Rabbit mAb
Remarks/Alias	p38; CRKII; CRK/CRKII
Species	Human
GeneID (human)	1398
GeneID	1398
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CRK/CRKII (P46108).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46108
Protein Weight	34kDa
Shipping	Ice bag