

Anti-Crkl (Phospho-Tyr207) Rabbit Polyclonal Antibody For Wb, Elisa - P46109

Item Number: CM0001054



Introduction

Phospho-CRKL-Y207 Rabbit pAb is a rabbit polyclonal antibody targeting phosphorylated CRKL (Tyr207) of human origin, validated for WB and ELISA with mouse cross-reactivity. Ideal for studying intracellular signal transduction.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CRKL-Y207 Rabbit pAb
Remarks/Alias	CRKL; Phospho-CRKL-Y207
Species	Human
GeneID	1399
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y207 of human CRKL (NP_005198.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P46109
Protein Weight	34kDa
Shipping	Ice bag