

Anti-Crkl Polyclonal Antibody For Wb, Elisa - P46109

Item Number: CM0014892



Introduction

Rabbit polyclonal antibody against human CRKL, validated for Western blot and ELISA. Reacts with human and mouse CRKL. Suitable for signal transduction research and diagnostic assay development.

[Learn More](#)

Parameter	Description
Product Name	CRKL Rabbit pAb
Remarks/Alias	CRKL
Species	Human
GeneID (Human)	1399
GeneID	1399
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46109
Protein Weight	34kDa
Shipping	Ice bag