

Anti-Crkl Rabbit Polyclonal Antibody For Wb, Elisa - P46109

Item Number: CM0017464



Introduction

CRKL rabbit polyclonal antibody validated for WB and ELISA. Detects Crk-like protein (34 kDa) in human, mouse, and rat samples. Suitable for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	CRKL Rabbit pAb
Remarks/Alias	CRKL
Species	Human
GeneID (Human)	1399
GeneID	1399
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-303 of human CRKL (NP_005198.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46109
Protein Weight	34kDa
Shipping	Ice bag