

Anti-Crlf2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hc73

Item Number: CM0015014



Introduction

Rabbit polyclonal antibody against CRLF2 (TSLP receptor), validated for WB and ELISA. Reacts with human, mouse, and rat samples. Suitable for hematopoietic system and cytokine signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CRLF2 Rabbit pAb
Remarks/Alias	CRL2; TSLPR; CRLF2Y; CRLF2
Species	Human
Gene ID (Human)	64109
Gene ID	64109
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-231 of human CRLF2 (NP_071431.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HC73
Protein Weight	42kDa
Shipping	Ice bag