

Anti-Crls1 Polyclonal Antibody For Wb, Elisa - Q9Uja2

Item Number: CM0015019

Introduction

CRLS1 Rabbit pAb validated for WB and ELISA against human cardiolipin synthase, a mitochondrial inner membrane enzyme crucial for cardiolipin synthesis. Recombinant immunogen, ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	CRLS1 Rabbit pAb
Remarks/Alias	CLS; CLS1; GCD10; COSPD57; C20orf155; dj967N21.6; CRLS1
Species	Human
Gene ID (Human)	54675
Gene ID	54675
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human CRLS1 (NP_061968.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UJA2
Protein Weight	33kDa
Shipping	Ice bag