


Anti-Crm1/Xpo1 Monoclonal Antibody For Wb, If/Icc, Elisa - O14980

Item Number: CM0003339



Introduction

Rabbit monoclonal antibody against human CRM1/XPO1 (Exportin-1). Suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying nuclear export and host-virus interactions. High-quality IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	CRM1/XPO1 Rabbit mAb
Remarks/Alias	emb; CRM1; exp1; CRM-1; Exportin 1; CRM1/XPO1
Species	Human
Gene ID (Human)	7514
Gene ID	7514
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 972-1071 of human CRM1/XPO1 (O14980).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14980
Protein Weight	123kDa
Shipping	Ice bag