

Anti-Crm1/Xpo1 Rabbit Polyclonal Antibody For Wb And Elisa - O14980

Item Number: CM0016232



Introduction

Rabbit polyclonal antibody targeting human Exportin-1 (CRM1/XPO1), validated for Western blot and ELISA, with mouse cross-reactivity. Suitable for nuclear export and viral-host interaction research.

[Learn More](#)

Parameter	Specification
Product Name	CRM1/XPO1 Rabbit pAb
Remarks/Alias	emb; CRM1; exp1; CRM-1; CRM1/XPO1
Species	Human
Gene ID (Human)	7514
Gene ID	7514
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14980
Protein Weight	123kDa
Shipping	Ice bag