


Anti-Crmp1 Rabbit Monoclonal Antibody For Wb And Elisa - Q14194

Item Number: CM0004780



Introduction

Rabbit monoclonal antibody targeting human CRMP1 (DRP1, CRMP-1), validated for WB and ELISA. Cross-reacts with rat; predicted MW ~62 kDa. Ideal for axon guidance and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	CRMP1 Rabbit mAb
Remarks/Alias	DRP1; DRP-1; CRMP-1; DPYSL1; ULIP-3; CRMP1
Species	Human
GeneID (Human)	1400
GeneID	1400
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 450-550 of human CRMP1 (Q14194).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q14194
Protein Weight	62 kDa
Shipping	Ice bag