


Anti-Crmp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14194

Item Number: CM0019045



Introduction

CRMP1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets 62kDa dihydropyrimidinase-related protein 1 involved in axon guidance and cytoskeleton remodeling.

[Learn More](#)

Parameter	Value
Product Name	CRMP1 Rabbit pAb
Remarks/Alias	DRP1, DRP-1, CRMP-1, DPYSL1, ULIP-3, CRMP1
Species	Human
Gene ID (Human)	1400
Gene ID	1400
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 460-540 of human CRMP1 (NP_001304.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14194
Protein Weight	62kDa
Shipping	Ice bag