


Anti-Crmp2/Dpysl2 Monoclonal Antibody For Wb, If-P, Elisa - Q16555

Item Number: CM0012200



Introduction

Rabbit monoclonal antibody against human CRMP2/DPYSL2 (UniProt Q16555). Validated for WB, IF-P, and ELISA. Reacts with human, mouse, and rat. Useful for neuronal development, axon guidance, and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CRMP2/DPYSL2 Rabbit mAb
Remarks/Alias	DRP2; N2A3; CRMP2; DRP-2; ULIP2; CRMP-2; DHPRP2; ULIP-2; L2
Species	Human
Gene ID (Human)	1808
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 473-572 of human CRMP2/DPYSL2 (Q16555)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16555
Protein Weight	62kDa
Shipping	Ice bag