


Anti-Crmp2/Dpysl2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16555

Item Number: CM0018289



Introduction

Rabbit polyclonal antibody targeting CRMP2/DPYSL2, validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Suitable for studying neuronal development, axon guidance, and CRMP2 function.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CRMP2/DPYSL2 Rabbit pAb
Remarks/Alias	DRP2; N2A3; CRMP2; DRP-2; ULIP2; CRMP-2; DHPRP2; ULIP-2; L2
Species	Human
Gene ID (Human)	1808
Gene ID	1808
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human CRMP2/DPYSL2 (NP_001377.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16555
Protein Weight	62kDa
Shipping	Ice bag