


Anti-Crmp3/Dpysl4 Monoclonal Antibody For Wb, Elisa - O14531

Item Number: CM0003414



Introduction

Recombinant rabbit monoclonal anti-CRMP3/DPYSL4 antibody validated for Western blot and ELISA. Recognizes human and rat species. Ideal for investigating axon guidance, neuronal growth cone collapse, and cell migration.

[Learn More](#)

Parameter	Value
Product Name	CRMP3/DPYSL4 Rabbit mAb
Remarks/Alias	CRMP3; DRP-4; ULIP4; CRMP3/DPYSL4
Species	Human
GeneID (Human)	10570
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CRMP3/DPYSL4 (O14531)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	O14531
Protein Weight	62kDa
Shipping	Ice bag