


Anti-Crmp5/Dpysl5 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Bpu6

Item Number: CM0010991



Introduction

Rabbit monoclonal antibody targeting CRMP5/DPYSL5 (Dihydropyrimidinase-related protein 5), validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for neuroscience studies on dendrite outgrowth regulation.

[Learn More](#)

Parameter	Value
Product Name	CRMP5/DPYSL5 Rabbit mAb
Remarks/Alias	CV2; CRAM; CRMP5; RTSC4; Ulip6; CRMP-5; CRMP5/DPYSL5
Species	Human
GeneID (human)	56896
GeneID	56896
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 465-564 of human CRMP5/DPYSL5 (Q9BPU6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BPU6
Protein Weight	61kDa
Shipping	Ice bag