

Anti-Crnn Rabbit Polyclonal Antibody (Wb, Elisa) - Q9Ubg3

Item Number: CM0030065

Introduction

CRNN Rabbit Polyclonal Antibody validated for Western blot and ELISA; reacts with human and mouse Cornulin. Ideal for studying cell proliferation and squamous epithelium biology. - Q9UBG3

[Learn More](#)

Parameter	Value
Product Name	CRNN Rabbit pAb
Remarks/Alias	DRC1; PDRC1; SEP53; C1orf10; CRNN
Species	Human
Gene ID (Human)	49860
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 396-495 of human CRNN (NP_057274.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UBG3
Protein Weight	54kDa
Shipping	Ice bag