


Anti-Crtc1 Rabbit Polyclonal Antibody For Wb And Elisa - Q6Uuv9

Item Number: CM0014375



Introduction

Anti-CRTC1 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human CRTC1 and cross-reacts with mouse and rat. Ideal for CREB signaling, L-LTP, and circadian clock research.

[Learn More](#)

Parameter	Value
Product Name	CRTC1 Rabbit pAb
Remarks/Alias	MAML2; MECT1; Mam-2; TORC1; TORC-1; WAMTP1; CRTC1
Species	Human
Gene ID (Human)	23373
Gene ID	23373
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CRTC1 (NP_056136.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6UUV9
Protein Weight	67kDa
Shipping	Ice bag