

Anti-Crtc2 Rabbit Polyclonal Antibody For Wb, Elisa - Q53Et0

Item Number: CM0014401



Introduction

Rabbit polyclonal antibody targeting human CRTC2 (TORC2), validated for Western blot and ELISA. Detects a 73kDa band; suitable for research on CREB signaling, gluconeogenesis, and metabolic regulation.

[Learn More](#)

Parameter	Specification
Product Name	CRTC2 Rabbit pAb
Remarks/Alias	TORC2; TORC-2; CRTC2
Species	Human
GeneID (human)	200186
GeneID	200186
Immunogen	Recombinant fusion protein containing amino acids 16-126 of human CRTC2 (NP_859066.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q53ET0
Protein Weight	73kDa
Shipping	Ice bag