

Anti-Crx Rabbit Polyclonal Antibody For Wb, If-P, Elisa - O43186

Item Number: CM0024761



Introduction

Rabbit polyclonal antibody against human Cone-rod homeobox protein (CRX). Validated for Western blot, immunofluorescence (paraffin), and ELISA. Cross-reacts with mouse. Ideal for photoreceptor gene regulation and retinal development research.

[Learn More](#)

Parameter	Value
Product Name	CRX Rabbit pAb
Remarks/Alias	CRD; LCA7; OTX3; CORD2; CRX
Species	Human
Gene ID (Human)	1406
Gene ID	1406
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CRX (NP_000545.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse
SWISS	O43186
Protein Weight	32kDa
Shipping	Ice bag