

Anti-Cry1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16526

Item Number: CM0013192

Introduction

Rabbit polyclonal antibody against human CRY1 (Cryptochrome-1), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for circadian rhythm studies. 66kDa, shipped on ice.

[Learn More](#)

Parameter	Value
Product Name	CRY1 Rabbit pAb
Remarks/Alias	DSPD; PHLL1; CRY1
Species	Human
GeneID (Human)	1407
GeneID	1407
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 496-586 of human CRY1 (NP_004066.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16526
Protein Weight	66kDa
Shipping	Ice bag