

Anti-Cry2 Rabbit Polyclonal Antibody For Wb, Elisa - Q49An0

Item Number: CM0029115

Introduction

High-quality rabbit polyclonal antibody targeting human CRY2, a key circadian clock repressor. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for circadian rhythm and metabolism research. Expected band at 67 kDa.

[Learn More](#)

Parameter	Value
Product Name	CRY2 Rabbit pAb
Remarks/Alias	HCRY2, PHLL2, CRY2
Species	Human
GeneID (Human)	1408
GeneID	1408
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human CRY2 (NP_066940.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q49AN0
Protein Weight	67kDa
Shipping	Ice bag