

Anti-Cry2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q49An0

Item Number: CM0015185

Introduction

Rabbit polyclonal antibody targeting human CRY2 (67kDa), validated for WB, IF/ICC, and ELISA. Recognizes human and rat CRY2. Suitable for circadian clock, transcription regulation, and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	CRY2 Rabbit pAb
Remarks/Alias	HCRY2; PHLL2; CRY2
Species	Human
GeneID (Human)	1408
GeneID	1408
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 455-614 of human CRY2 (NP_066940.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q49AN0
Protein Weight	67kDa
Shipping	Ice bag