


Anti-Cryaa Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P02489

Item Number: CM0010680



Introduction

Rabbit monoclonal antibody against human CRYAA (alpha-crystallin A chain/HspB4). Validated for WB, IF-P, and ELISA; cross-reacts with mouse and rat. Ideal for lens biology and chaperone function studies.

[Learn More](#)

Parameter	Value
Product Name	CRYAA Rabbit mAb
Remarks/Alias	CRYA1; HSPB4; CTRCT9
Species	Human
GeneID (Human)	1409
GeneID	1409
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 6-95 of human CRYAA (P02489)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P02489
Protein Weight	20kDa
Shipping	Ice bag