


Anti-Cryaa Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P02489

Item Number: CM0028545



Introduction

Rabbit polyclonal antibody targeting Alpha-crystallin A chain (CRYAA; HSPB4) validated for WB, IF-P, and ELISA in human, mouse, rat, zebrafish. Suitable for lens biology and chaperone function research.

[Learn More](#)

Parameter	Value
Product Name	CRYAA Rabbit pAb
Remarks/Alias	CRYA1; HSPB4; CTRCT9; CRYAA
Species	Human
Gene ID	1409 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-173 of human CRYAA (NP_000385.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Zebrafish
SWISS	P02489
Protein Weight	20kDa
Shipping	Ice bag