


Anti-Cryab Polyclonal Antibody For Wb And Elisa - P02511

Item Number: CM0018068



Introduction

Anti-CRYAB polyclonal antibody (Rabbit) validated for WB and ELISA, reacting with human, mouse, and rat CRYAB. Ideal for research on crystallin biology and stress responses.

[Learn More](#)

Parameter	Value
Product Name	CRYAB Rabbit pAb
Remarks/Alias	MFM2; CRYA2; CTPP2; HSPB5; CMD1II; CTRCT16; HEL-S-101; CRYAB
Species	Human
GeneID (Human)	1410
GeneID	1410
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human CRYAB (NP_001876.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02511
Protein Weight	20kDa
Shipping	Ice bag